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
NOTICE TO DISCONNECT THE TRANSFER PART

The transfer part of OKTAGON® implants is used to transfer force from the insertion tool to the implant during insertion. For constructional reason it is not possible to completely exclude the possibility of jamming between the implant and the transfer part, especially at higher torques than the recommended 35 Ncm. This means that the transfer part might not simply be removed vertically from the implant. We would like to suggest the following course of action.

In case of jamming between implant and transfer part, please follow the steps below:

1. **Manual loosening of the transfer part using a hand ratchet**
 - a. Limit the insertion torque of the hand ratchet to 10 Ncm (or $\frac{1}{2}$ x insertion torque)
 - b. Plug OKTAGON® adapter for ratchet (ref. 64558 or 17372) onto transfer part
 - c. Place hand ratchet in „OUT“ setting on the insertion tool
 - d. Carefully turn counterclockwise until a maximum click is heard from the hand ratchet.
The transfer part can then be removed from the implant by pulling it off vertically
2. **Mechanical loosening of the transfer part by means of the drive unit**
 - a. Set the insertion torque limit on the drive unit to 10 Ncm
(or $\frac{1}{2}$ x insertion torque)
 - b. Set the direction of rotation of the handpiece to counterclockwise
 - c. Carefully turn back and release the transfer part vertically from the implant

NOTE: When releasing the clamping according to the points described above, only the transfer part is turned back by the torque limiter, not the implant.

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